
	Document Number:	SOP-101-A_2
	Effective Date:	08/24/2022
	Version:	2
	1X Inactivation Saline Solution Prep - Table Form	

<u>100X Inactivation Solution</u> Lot (-MMDD): Freeze EXP (MM-YY): Thaw EXP (+3 months): <u>Put New in Use:</u> Location: Name: Date and Time:	<div style="border: 1px solid black; padding: 10px; text-align: center;"> QR CODE </div> Log Link: SOP-101-B_1 {SITE} Name: Date and Time:
---	--

1X Inactivation Saline Solution Prep Protocol:

- 1) Vortex 100X Inactivation Solution for 10 sec.
- 2) Add 1/100 volume of 100X Inactivation Solution to Saline.
- 3) Vortex for 10 sec, Invert 3 times, and Vortex for 10 more sec.
 * If making Inactivation Solution in a bottle, shake instead of vortexing
- 4) Label tube with 1XISS ID (1X_MMDD-1).
- 5) EXPIRATION of 1X Inactivation Saline Solution: 2 weeks after prep.

Name	Date + Time	1XISS ID (1X_MMDD-1)	Saline Lot + Exp	Saline mL	100XISS µL

	Document Number:	SOP-101-A_2
	Effective Date:	08/24/2022
	Version:	2
	1X Inactivation Saline Solution Prep - Table Form	

Name	Date + Time	1XISS ID (1X_MMDD-1)	Saline Lot + Exp	Saline mL	100XISS µL