



Collection Kit Assembly Instructions

Before creating kits

- Clean all work surfaces with bleach (clorox wipes, etc.) and ethanol/isopropyl alcohol
- Wash hands thoroughly

While creating kits

- Wear a mask (N95 ideal)
- Wear gloves (make sure not to touch gloves's fingers with bare hands)

Kit contents

- 1 slider top bag with the following inside:
 - 4 tubes
 - 16 swabs
 - 4 bio-bags
 - Instruction pamphlet
 - Sticker labels with family names
 - Place 1 sticker on outside of slider top bag
 - Place the remaining stickers inside the bag facing out

